7921RA: BATH SEAT WITH ARMS AND BACKREST / 7931RA: BATH SEAT WITH ARMS

ASSEMBLY AND OPERATION INSTRUCTIONS

SAVE THESE INSTRUCTIONS FOR FUTURE USE

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IMPORTANT SAFETY INFORMATION: PLEASE READ

- ⚠ **WARNING:** Important! Read and understand these instructions before assembling or using the Bath Seat. If you do not understand any part of these warnings, cautions or instructions, contact a healthcare professional for direction in the use of this product. If the Bath Seat is not properly assembled and adjusted, personal injury and damage to the Bath Seat could result.
- ⚠ **WARNING:** If components are damaged or missing, contact your distributor immediately. **DO NOT** use substitute parts. Use only Lumex replacement parts. Non-Lumex replacement parts could cause personal injury and damage to the Bath Seat.
- ⚠ **WARNING:** Cancer and Reproductive Harm - www.p65warnings.ca.gov.
- ⚠ **WARNING:** GF Health Products, Inc. assumes no responsibility for any damage or injury caused by improper assembly or use of this product.


INTENDED USE

The 7921RA / 7931RA Bath Seat is intended for use by individuals who require seating support while bathing or showering. It is designed to be placed fully inside the bathtub or shower.

Contraindications

- ⚠ **WARNING:** This product has a weight limit of 350 lb, **EVENLY DISTRIBUTED**.
- ⚠ **WARNING:** **DO NOT** use this product as a transfer bench or climbing device.
- ⚠ **WARNING:** Use this product only with tub floors wider than 18 3/4".

ASSEMBLY

7921RA / 7931RA Bath Seat with Arms
(7921RA with Backrest Shown Above)

1. Remove all components (shown above) from the box.

2. Place the seat facedown on a flat surface.
3. Locate the four legs.
4. Install legs on frame: press lock button on each plastic leg sleeve end and slide the leg onto the frame until the button clicks into place.
5. Adjust leg extension height by pushing in the leg height adjustment lock button, moving to the desired hole height, and releasing button until it clicks into place.


Info: *If assembling 7921RA (with backrest), proceed to step 6.*

If assembling 7931RA (without backrest), proceed to step 9.

6. Locate backrest tubes and backrest.
7. Install backrest tubes on backrest: press each backrest tube lock button and slide tube onto backrest until button clicks into place.
8. Install backrest on bath seat: press both backrest lock buttons and slide backrest tubes into bath seat backrest channel guides until buttons click into place.
9. Attach arms: press both arm lock buttons and insert each arm into arm socket until arm clicks into place.

⚠ WARNING: Ensure all lock buttons have clicked into place and all lock buttons protrude fully through holes. Contact your GF authorized distributor immediately if you see any issue with fit, form, function, or safety.

⚠ WARNING: Ensure the Bath Seat is assembled as described. If the Bath Seat is not properly assembled and installed, personal injury and damage to the Bath Seat could result.

**7921RA Bath Seat
with Arms and Backrest**

OPERATION

Adjusting Leg Height

Push in the leg height adjustment lock button and move leg extension to the desired height position — you will hear it click into place. Ensure that all four legs are locked into place at the same height and all leg height adjustment lock buttons protrude fully through adjustment holes.

⚠ WARNING: Ensure all Bath Seat components are secure before each use.

Center of Balance During Use

The point where the Bath Seat will tip forward, rearward or sideward is determined by its center of balance and stability, and the center of balance and stability is affected by the user's body position and weight distribution. If you fail to heed the following warnings, falling or tipping may occur and could cause severe personal injury.

⚠ WARNING: To avoid falling or tipping:
Have someone help you until you know the balance points of your Bath Seat and how to avoid tipping.

Do not put downward weight on one armrest at a time. If using the armrest for balance, hold both armrests at once to distribute your weight equally.

Do not put downward weight on the backrest.

Do not reach or lean if you must shift your weight to do so. Ask for help or use a device to extend your reach.

Never try to reach an object if you must move forward in your seat and your buttocks do not stay in full contact with the seat.

Never lean backward to reach an object.

MAINTENANCE

⚠ WARNING: Check the Bath Seat weekly for signs of damage or cracking. If damage or cracking of any component of the Bath Seat is found, DO NOT use the Bath Seat. Contact your Graham-Field authorized provider.

To clean the Bath Seat, wash it with soap and water or a mild detergent, then rinse. Ensure the Bath Seat is rinsed well before use.

**7921RA Bath Seat
with Arms and Backrest
Shown Assembled**

**7931RA Bath Seat with Arms
Shown Assembled**

LIMITED WARRANTY

SCOPE OF WARRANTY

GF Health Products, Inc. ("GF") warrants to the original purchaser only that it will replace or repair components, at GF's sole discretion, that are defective in material or workmanship under normal use and service. All warranties are conditioned upon the proper use of the products strictly in accordance with good commercial practice and applicable GF instructions and manuals, including proper use and maintenance. To the extent that a component is warranted by a third party, GF conveys all of its rights under that warranty to the original purchaser, to the extent permitted. This limited warranty shall only apply to defects that are reported to GF's customer service team within the applicable warranty period and which, upon examination by GF or its authorized representative, prove to be a warranty item. This limited warranty is not transferable.

The warranted components and time period are set forth below:

- Frame.....three (3) years
- Seat/back, leg extensions, and arms.....six (6) months
- Normal wear components such as tips, wheels/casters, and handgrips.....three (3) months

The applicable warranty period shall commence from date of shipment to the original customer, unless there is an expiration date on the component in which case the warranty shall expire on the earlier of warranty period or the expiration date.

OBTAINING WARRANTY SERVICE

This limited warranty shall only apply to defects that are reported to the Distributor from whom the Customer purchased the product within the applicable warranty period. If there is not a Distributor, you must contact GF directly by calling 1.770.368.4700, sending a fax request to 1.770.368.2386, or by e-mailing a request to cs@grahamfield.com. Specific directions will be provided by the Customer Service Representative. Failure to abide by the specific directions will result in denial of the warranty claim.

EXCLUSIONS

The warranty does not cover and GF shall not be liable for the following:

- 1) Defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- 2) Products which are not installed, used, or properly cleaned and maintained as required in the official manual for the applicable product;
- 3) Products considered to be of a non-durable nature including, but not limited to: filters, fuses, gaskets, lubricants, and charts;
- 4) Accessories or parts not provided by GF;
- 5) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing, in advance, by GF;
- 6) Any labor or shipping charges incurred in the replacement part installation or repair;
- 7) Costs and expenses of regular maintenance and cleaning; and
- 8) Representations and warranties made by any person or entity other than GF.

ENTIRE WARRANTY, EXCLUSIVE REMEDY AND CONSEQUENTIAL DAMAGES DISCLAIMER

THIS WARRANTY IS GF'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IF ANY MODEL OR SAMPLE WAS SHOWN TO THE CUSTOMER, SUCH MODEL OR SAMPLE WAS USED MERELY TO ILLUSTRATE THE GENERAL TYPE AND QUALITY OF THE PRODUCT AND NOT TO REPRESENT THAT THE PRODUCT WOULD NECESSARILY CONFORM TO THE MODEL OR SAMPLE IN ALL RESPECTS. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF THE DEFECTIVE PARTS. GF SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO: DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER, OR EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS.

The warranties contained herein contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document.

For additional information on this product or this warranty, please contact a GF Customer Service Representative.

NOTES:

- 1) Additional terms and conditions may apply.
- 2) Freight claims must be notated on the appropriate shipping documents and must be made with immediacy. International, federal and state regulations govern specific requirements for freight claims. Failure to abide by those regulations may result in a denial of the freight claim. GF will assist you in filing the freight claim.
- 3) Claims for any short shipment must be made within three (3) days of the invoice date.

1.770.368.4700

Information contained herein is subject to change.
The most current and complete product information can be found on our website.
www.grahamfield.com

© 2010, GF Health Products, Inc. All Rights Reserved.
Lumex is a trademark of GF Health Products, Inc.,
One Graham-Field Way, Atlanta GA 30340-3140
GF Health Products, Inc. is an ISO 13485 Certified Company.
Manufactured for GF Health Products, Inc.
Made in China

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© 2010, GF Health Products, Inc. Todos los derechos reservados.
Lumex es una marca de GF Health Products, Inc.,
GF Health Products, Inc. es una empresa certificada ISO 13485.
Fabricado para GF Health Products, Inc.