

MOBILITY ACCESSORIES

GLIDE BRAKES

Features

- For use in conjunction with swivel wheels or fixed wheels, glide brakes provide a safe, gradual braking action for rear walker legs
- Glide brakes attach to the rear of walker legs
- By pressing down on the handgrip, glide brakes allow rear walker tips to gradually touch the surface preventing forward motion
- Designed to meet the requirements of HCPCS code: E0159
- Limited lifetime warranty

Glide Brakes are recommended for use on Lumex Walkers. (We cannot guarantee the performance of these accessories on other industry products.)



		UOM
603910A	Glide Brakes	1 pr