



MOBILITY ACCESSORIES

PLATFORM WALKER ATTACHMENT

Features

- Designed for users with limited strength or dexterity who are unable to grip the walker
- Adjustable strap comfortably secures arm in the forearm rest
- Padded forearm rest can be pivoted and angled to fit user comfortably and to assist in lifting walker
- Padded handgrip provides comfortUniversal clamp fits most walkers
- Tool-free installation and adjustment allows quick fitting for most patient sizes and types
- Designed to meet the requirements of HČPCS code: E0154
- Three (3) year limited warranty





Platform Attachment is recommended for use on Lumex Walkers. (We cannot guarantee the performance of this product on other industry walkers.)



		UOM
6132B	Platform Walker Attachment	1 pr

* It is recommended that the platform walker attachments be installed on both sides for better stability.





www.grahamfield.com

